

# A New Lipid Software Workflow for Processing Orbitrap-based Global Lipidomics Data in Translational and Systems Biology Research

David A. Peake<sup>1</sup>, Yasuto Yokoi<sup>2</sup>, Junhua Wang<sup>1</sup> and Yingying Huang<sup>1</sup>

<sup>1</sup>Thermo Fisher Scientific, San Jose, USA; <sup>2</sup>Mitsui Knowledge Industry, Tokyo, JAPAN

## Overview

**Purpose:** We present a new workflow for high-resolution Thermo Scientific™ Orbitrap™ based mass spectrometers for lipidomics using a model system consisting of a wild-type strain vs. knockout for Co-Q production in yeast<sup>1</sup>.

**Methods:** Lipids in yeast mitochondria were analyzed by high resolution LC-MS and MS/MS. Lipid Search® software, an MS<sup>2</sup> based search using a comprehensive lipid database, was used to identify the lipid species and determine significant differences.

**Results:** The yeast lipidomics results obtained from the LC-MS data using Lipid Search are comparable to results obtained using infusion lipidomics. We also compared the lipids identified using metabolomics analysis of the same data set – component finding and molecular weight (MW) search for assignment of metabolites and lipids. Due to the complexity of lipid extracts we found that the comprehensive lipid database MS<sup>2</sup> search method is superior to the accurate mass based MW search for lipidomics.

## Introduction

Application of lipidomics to disease phenotype analysis is a growing area in research. Identification of unique biomarkers to distinguish healthy humans compared to individuals with disease can have an impact on the early detection of diseases and personalized medicine.

The complexity of the lipidome (Table 1) includes 8 major categories of lipids, over 80 major classes, 300 sub-classes and thousands of lipid species<sup>2</sup> many with overlapping isomeric or isobaric molecular ions. Because of this complexity, MW searches alone are not sufficient to identify lipids in a complex biological extract.

Identification of lipids requires sophisticated software with an extensive database. The combination of ultra-high resolution MS and MS<sup>n</sup> analysis should provide unambiguous and precise identification of lipids in biological samples. A robust algorithm for database searching of high-resolution data was developed by Prof. Ryo Taguchi and co-workers<sup>3</sup> and was commercialized by MKI (Tokyo, JP) as described recently<sup>4</sup>.

## Method

### Phenotypes of wild-type (WT) and Knockout (KO) Yeast Strains (*S. Cerevisiae*)

WT yeast continue to grow after glucose is exhausted from the media (Diauxic shift point, Figure 1) whereas KO yeast have a defect in Coenzyme Q production and do not grow after the shift. Duplicate biological replicates of WT and KO yeast were collected post shift (green triangle) for metabolomic/lipidomic analyses and analyzed by LC-MS.

### Sample Preparation

Yeast were treated with zymolase, homogenized and mitochondria were enriched by differential centrifugation. Mitochondrial protein levels were determined by BCA assay. Mitochondria (~0.25 mg) were extracted 3 times with 400 µL of IPA for 10 min at 4°C. After centrifugation, supernatants were combined and vacuum dried. Samples were dissolved in 250 µL of 65:35:5 Acetonitrile, Isopropanol, Water with 5 µg/mL 17:0 PG.

### Liquid Chromatography–Mass Spectrometry (LC-MS)

Thermo Scientific™ Accela™ 1250 chromatograph, Accela Open auto-sampler, 10 µL Injection. Column: 2.1 x 100 mm C18, 2.7µm, 55°C, 260 µL/min. HPLC method<sup>1</sup> is described in S. Bird, et al., *Anal. Chem.* **2011**, *83*, 940–949, 6648–6657. A Thermo Scientific™ Q Exactive™ high-resolution Orbitrap mass spectrometer was operated at 70K resolution for ESI pos. ion LC-MS and 35K for Top5 MS/MS (CE 35).

### Data Analysis Software

Metabolomics –Thermo Scientific™ SIEVE™ and Lipidomics – Thermo Scientific™ Lipid Search™.


**Table 1. Lipid Complexity from the LIPID MAPS Structure Database (LMSD)<sup>2</sup>**

Lipid Category	# Class	# Sub-Class	# Lipids
FA	14	36	5,787
GL	6	19	7,568
GP	21	120	8,001
SP	10	31	4,317
ST	6	38	2,678
PR	5	21	1,200
SL	6	7	1,293
PK	15	28	6,741
<b>Total</b>	<b>83</b>	<b>300</b>	<b>37,585</b>


## Results

### High-Resolution LC-MS Data – Metabolomics Analysis


To characterize the yeast phenotypes we analyzed the sample extracts using an LC-MS method suited for analysis of both metabolites and lipids. The LC-MS chromatogram from WT yeast (Figure 2) shows the regions where lipid classes elute during the LC gradient. Metabolomics analysis using an accurate-mass search tentatively identified 160 metabolites and lipids were present. Principal component analysis and t-Test statistics (Figure 3) showed many key metabolite differences.



**Figure 1. Growth Phenotypes of WT and KO Yeast Strains.**




**Figure 2. LC-MS Chromatograms of Lipids from WT and KO Yeast.**



**Figure 3. Significant Metabolite Differences Observed WT vs. KO Yeast.**

### LC-MS/MS Data Processing Workflow using Lipid Search Software (Figure 4)

- 1. Peak Detection.** Read raw files, MS<sup>n</sup> and precursor ion accurate masses.
- 2. Identification.** Candidate molecular species are identified by searching a large database > 1,500,000 entries of accurate masses (lipid precursor and fragment ions) predicted from each potential lipid structure and positive / negative ion adducts.
- 3. Alignment.** The search results for each individual sample are aligned within a time window and the results are combined into a single report.
- 4. Quantification.** The accurate-mass extracted ion chromatograms are integrated for each identified lipid precursor and the peak areas are obtained.
- 5. Statistical Analysis.** t-Tests determine which lipid species are significantly different between sample vs. control groups, and results are displayed in a whisker plot.



**FIGURE 4. Lipid Search Software LC-MS Workflow.**

### Submitting Data for Lipid Search Identification and Alignment

LC-MS raw data files containing full scan and data dependent-MS/MS were searched for PL, GL, SP and Co-enzyme lipid classes using a mass tolerance of 5 ppm for precursor ions and 8 ppm for product ions (Figure 5a).

The search results from the 4 samples were aligned using a 0.25 min tolerance window and a combined report was generated (Figure 5b).

Name	RawData	Type	ExpType	Process	Result	Request	Update
11WTPost_1.raw	Product	LC	I	Q	322	2013/10/08 11:00:15	2013/10/08 11:07:27
12WTPost_1.raw	Product	LC	I	Q	322	2013/10/08 11:00:15	2013/10/08 11:06:55
3KOPost_1.raw	Product	LC	I	Q	309	2013/10/08 11:00:15	2013/10/08 11:07:22
4KOPost_1.raw	Product	LC	I	Q	306	2013/10/08 11:00:15	2013/10/08 11:06:52

**FIGURE 5a. Search Results for Yeast Lipids.**


Name	Type	ExpType	Process	Result	Request	Update
Yeast (All) Q	Product	LC	M	738	2013/10/10 09:47:18	2013/10/10 09:48:54
Yeast (All) Q	Product	LC	M	542	2013/10/10 09:55:45	2013/10/10 09:57:08

**FIGURE 5b. Alignment Results for Yeast Lipids.**

Search results obtained in < 10 min with 64-bit laptop (MS Windows 7, 2.2 GHz, Intel i7 CPU, 8GB RAM).

### Identification Report (Figure 6)


For each MS<sup>2</sup> spectrum, search results are summarized for lipid species matching the predicted fragmentation pattern from the database with a score indicating the fit. If a mixture of lipids is found, the most abundant lipid is displayed. The fragment ions used to identify the lipid are highlighted in red when each of the species are selected.



**Figure 6. Search Results for m/z 584.5249, Rt = 17.3 min, DG (32:1).**

### Combined Report – Details (Figures 7 and 8)

Lipid species identified in each LC-MS data file were aligned across the dataset within a retention time tolerance. Quantification is performed on the relative amount of the precursor ion, which in some cases was identified as a mixture of isomers. For each lipid species in the aligned dataset, an interactive report allows review of the data. Relative amounts of each identified lipid were quantified by peak areas and significant differences were determined using t-Tests (Table 2) producing a heat map.



**FIGURE 7. Combined Report Results for PG (17:0/17:0) Internal Standard.**



**Figure 8. Combined Report Results – Total Lipid Profile.**

### Yeast Lipidomics Results

The total number of sum composition lipids identified in yeast WT and KO mitochondria (542) is twice the number of lipids quantified previously (250) by infusion lipidomics<sup>5</sup>.

**Table 2. Summary of Differences between WT vs. KO Yeast Lipids. Analyses with p-Values < 0.05 for t-Test between WT and KO groups. Fold-change (KO vs. WT) indicated by Red (increase) or Green (decrease).**

Class	Compound	RT min	Ratio	p-Value	Class	Compound	RT min	Ratio	p-Value	Class	Compound	RT min	Ratio	p-Value
Cer	Cer(d18:0/18:0)	16.73	0.006	0.007	Cer	Cer(d16:0/16:0)	15.45	0.007	0.027	DG	DG(16:0/16:0)	16.54	0.007	0.027
	Cer(d18:0/18:0)	15.06	0.006	0.007		Cer(d16:0/16:0)	16.54	0.007	0.027		DG(16:0/16:0)	16.54	0.007	0.027
	Cer(d18:0/18:0)	18.77	0.037	0.008		Cer(d16:0/16:0)	14.84	0.007	0.009		DG(16:0/16:0)	14.81	0.007	0.048
SP	SP(d18:0/18:0)	22.74	0.011	0.028	SP	SP(d18:0/18:0)	20.53	0.004	0.004	PE	PE(16:0/16:0)	19.54	0.009	0.009
	SP(d18:0/18:0)	22.69	0.002	0.004		SP(d18:0/18:0)	22.65	0.002	0.002		PE(16:0/16:0)	22.65	0.002	0.002
	SP(d18:0/18:0)	16.72	0.049	0.002		SP(d18:0/18:0)	20.01	0.002	0.002		PE(16:0/16:0)	20.01	0.002	0.002
GP	GP(d18:0/18:0)	3.03	0.026	0.002	GP	GP(d18:0/18:0)	32.25	0.002	0.002	TG	TG(16:0/16:0)	23.25	0.002	0.002
	GP(d18:0/18:0)	4.79	0.003	0.003		GP(d18:0/18:0)	32.55	0.002	0.002		TG(16:0/16:0)	23.25	0.002	0.002
	GP(d18:0/18:0)	10.02	0.000	0.000		GP(d18:0/18:0)	32.55	0.002	0.002		TG(16:0/16:0)	23.25	0.002	0.002
PE	PE(16:0/16:0)	22.62	0.004	0.004	PE	PE(16:0/16:0)	13.49	0.008	0.008	TG	TG(16:0/16:0)	23.25	0.002	0.002
	PE(16:0/16:0)	12.88	0.000	0.000		PE(16:0/16:0)	13.49	0.008	0.008		TG(16:0/16:0)	23.25	0.002	0.002
	PE(16:0/16:0)	14.36	0.045	0.004		PE(16:0/16:0)	13.89	0.009	0.009		TG(16:0/16:0)	23.25	0.002	0.002
GP	GP(16:0/16:0)	9.58	0.029	0.002	GP	GP(16:0/16:0)	14.26	0.003	0.003	TG	TG(16:0/16:0)	22.28	0.002	0.002
	GP(16:0/16:0)	11.78	0.007	0.007		GP(16:0/16:0)	14.26	0.003	0.003		TG(16:0/16:0)	22.28	0.002	0.002
	GP(16:0/16:0)	12.29	0.024	0.004		GP(16:0/16:0)	14.26	0.003	0.003		TG(16:0/16:0)	22.28	0.002	0.002
GP	GP(16:0/16:0)	15.71	0.029	0.002	GP	GP(16:0/16:0)	14.26	0.003	0.003	TG	TG(16:0/16:0)	22.28	0.002	0.002
	GP(16:0/16:0)	12.70	0.015	0.002		GP(16:0/16:0)	14.26	0.003	0.003		TG(16:0/16:0)	22.28	0.002	0.002
	GP(16:0/16:0)	18.37	0.021	0.002		GP(16:0/16:0)	14.26	0.003	0.003		TG(16:0/16:0)	22.28	0.002	0.002
GP	GP(16:0/16:0)	12.40	0.000	0.000	GP	GP(16:0/16:0)	14.26	0.003	0.003	TG	TG(16:0/16:0)	22.28	0.002	0.002
	GP(16:0/16:0)	12.70	0.000	0.000		GP(16:0/16:0)	14.26	0.003	0.003		TG(16:0/16:0)	22.28	0.002	0.002
	GP(16:0/16:0)	12.40	0.000	0.000		GP(16:0/16:0)	14.26	0.003	0.003		TG(16:0/16:0)	22.28	0.002	0.002
GP	GP(16:0/16:0)	12.70	0.000	0.000	GP	GP(16:0/16:0)	14.26	0.003	0.003	TG	TG(16:0/16:0)	22.28	0.002	0.002
	GP(16:0/16:0)	12.70	0.000	0.000		GP(16:0/16:0)	14.26	0.003	0.003		TG(16:0/16:0)	22.28	0.002	0.002
	GP(16:0/16:0)	12.40	0.000	0.000		GP(16:0/16:0)	14.26	0.003	0.003		TG(16:0/16:0)	22.28	0.002	0.002
GP	GP(16:0/16:0)	11.85	0.000	0.000	GP	GP(16:0/16:0)	13.48	0.000	0.000	TG	TG(16:0/16:0)	22.28	0.002	0.002
	GP(16:0/16:0)	14.24	0.002	0.002		GP(16:0/16:0)	13.97	0.005	0.007		TG(16:0/16:0)	22.28	0.002	0.002
	GP(16:0/16:0)	12.71	0.007	0.007		GP(16:0/16:0)	13.97	0.005	0.007		TG(16:0/16:0)	22.28	0.002	0.002
GP	GP(16:0/16:0)	11.85	0.000	0.000	GP	GP(16:0/16:0)	13.97	0.005	0.007	TG	TG(16:0/16:0)	22.28	0.002	0.002
	GP(16:0/16:0)	14.24	0.002	0.002		GP(16:0/16:0)	13.97	0.005	0.007		TG(16:0/16:0)	22.28	0.002	0.002
	GP(16:0/16:0)	12.71	0.007	0.007		GP(16:0/16:0)	13.97	0.005	0.007		TG(16:0/16:0)	22.28	0.002	0.002
GP	GP(16:0/16:0)	11.85	0.000	0.000	GP	GP(16:0/16:0)	13.97	0.005	0.007	TG	TG(16:0/16:0)	22.28	0.002	0.002
	GP(16:0/16:0)	14.24	0.002	0.002		GP(16:0/16:0)	13.97	0.005	0.007		TG(16:0/16:0)	22.28	0.002	0.002
	GP(16:0/16:0)	12.71	0.007	0.007		GP(16:0/16:0)	13.97	0.005	0.007		TG(16:0/16:0)	22.28	0.002	0.002
GP	GP(16:0/16:0)	11.85	0.000	0.000	GP	GP(16:0/16:0)	13.97	0.005	0.007	TG	TG(16:0/16:0)	22.28	0.002	0.002
	GP(16:0/16:0)	14.24	0.002	0.002		GP(16:0/16:0)	13.97	0.005	0.007		TG(16:0/16:0)	22.28	0.002	0.002
	GP(16:0/16:0)	12.71	0.007	0.007		GP(16:0/16:0)	13.97	0.005	0.007		TG(16:0/16:0)	22.28	0.002	0.002
GP	GP(16:0/16:0)	11.85	0.000	0.000	GP	GP(16:0/16:0)	13.97	0.005	0.007	TG	TG(16:0/16:0)	22.28	0.002	0.002
	GP(16:0/16:0)	14.24	0.002	0.002		GP(16:0/16:0)	13.97	0.005	0.007		TG(16:0/16:0)	22.28	0.002	0.002
	GP(16:0/16:0)	12.71	0.007	0.007		GP(16:0/16:0)	13.97	0.005	0.007</					